

### ZT-11 (Gun Type)

Non Contact Infrared Thermometer



Measures Forehead and Objects



Fast Measurement: less than 2 seconds



Large LCD | Memory: 50



Measuring Distance: 1-5cm



°C and °F Alternatives and Auto Power Off



Fever Warning with Red back light over 38.5°C/101.3°F

## SPECIFICATIONS

### Product Information

- Monitor Size: 86\*39\*136 mm (3.38\*1.54\*5.35 inches)
- LCD Size: 26\*23mm(1.02\*0.90 inches)
- Measuring Range:  
Body mode 34°-43°(93.2-109.4°F);  
Object mode 10°-50.5°(50-122.0°F)
- Display Resolution for body mode  
±0.2°C (35°-42°), ±0.4°F (95.0-107.6°F)

### Packaging Information

- 1 pc per box
- Box Size: 100\*48\*155 mm (4\*1.9\*6.1 inches)
- Manual included
- 2 AAA Batteries

\*Product Models, Designs and Specifications may vary from time to time. The same will be explicitly confirmed at the time of each order

Zovec Healthcare Private Limited | 303E, Sector-59, Faridabad - 121004, Haryana, INDIA | Ph: +91 129 4300777