



ZN-C1

Compressor Type Nebulizer



Excellent for treatment of Asthma, Bronchitis & other respiratory disorders



-  Long Life & Durable
-  Excellent Vibration Reduction
-  Strong Compressor with efficient performance

SPECIFICATIONS

Power	220V/50HZ	Air Flow	> 6l/m
Power Consumption	170 VA Max	Operating Temperature	10° C to 40° C
Medication Capacity	8ml	Operating Humidity	< 85% RH
Partical Size	0.5 to 5µm	Storage Temperature	-10°C to 70° C
MMAD	2.5µm	Storage Humidity	<85% RH
Noise Level	60 dBA at 1 meter	Dimension (L*W*H)	195 X 150 X 100 mm
Nebulisation Rate	> 0.35ml/min	Weight	1.5kg
Operating Pressure Range	0.086 to 1.106 Mpa		
Compressor Pressure Range	> 0.2 Mpa		

Nebulizer is a drug delivery device used to administer medication in the form of a mist inhaled into the lungs.

*Product Models, Designs and Specifications may vary from time to time. The same will be explicitly confirmed at the time of each order