



## ZB SERIES

Hot Water Bottle

- Beneficial for sore backs stomach aches and muscular tension.
- Alleviate stress

# ZB-HW1

## SPECIFICATIONS

### Key features

- Zovec hot water bottle made with high quality odourless materials.
- Designed with high quality rubber material.
- Easy to clean, colour-ful, leakproof material provide relief from aches in joints and pain associated with rheumatism.
- It is made of natural rubber. It retains heat for longer time period. Ribbed outer surface helps to maintain the liquid temperature.

### Technical Specifications

- Weight- 330 gms
- Thickness-1.6 mm
- Rubber Content - 56% (High Concentration of real rubber content - 56%.) One mold manufacturing makes bottle soft and smooth to touch, it stays warm for longer time as rubber is great insulator of heat.
- BPA and phthalate free material - its eco friendly and UL tested for safety and reliability.

### OTHER ZOVEC PRODUCTS



Pulse Oximeter



BP Monitor



Bathroom Scales



Thermometer



Nebulizer

Zovec Healthcare Private Limited

303E, Sector-59, Faridabad-121004, INDIA | +91 0129 4300777 | Info@zovechealthcare.com | www.zovechealthcare.com