



# ZT SERIES

Digital Thermometer



Thermometer Color Range



Display Box

## ZT-D1



Accurate Memory



Oral  
Axilla Rectal



Fast



Waterproof



Button Batteries



Buzzer Remind



Auto off



Storage Box

## SPECIFICATIONS

Operating Temperature Condition	10°C ~ 40°C 50°F ~ 104°F
Measurement Range	89.6°F ~ 104°F
Accuracy	± 0.2°F
Battery	AG3 or LR41
Battery Life	2 Years
Unit Weight	Approx 9.8g including batteries

Dimensions	125 X 19 X 11 mm
Beep Sound Functionality	10 long beeps in 10 secs for normal temperature 30 short beeps in 10 secs in case of high temperature
Auto-Power Off	After 8 minutes

\*Product Models, Designs and Specifications may vary from time to time. The same will be explicitly confirmed at the time of each order